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*All the information given in this manual is what available at the time of printing.
The manufacturer reserves the right to make changes to the product at any time without notice.
It is advisable to verify if there are any updates.*

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All rights reserved.

*Forbidden **any form** of print, duplication, reproduction and publication of this instruction manual, or
part of it, without the written agreement of GIBERTINI ELETTRONICA Srl*

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1.	TECHNICAL CHARACTERISTICS SUPER ALCOMAT	22
2.	INSTALLATION OF SUPER ALCOMAT	23
2.1	INSTALLATION PROCEDURE	23
3.	DESCRIPTION OF SUPER ALCOMAT	24
3.1	COMMAND PANEL	24
4.	IGNITION OF SUPER ALCOMAT	25
5.	MEASUREMENTS WITH SUPER ALCOMAT	26
5.1	RELATIVE DENSITY AT 20°C OF A DISTILLATE AND ASV (ALCOHOLIC STRENGTH BY VOLUME)	27
5.1.1	Relative density of a liquid	28
5.1.2	Correction of the floater or compensation	28
5.1.2.1	Compensation of the floater with the hydroalcoholic mixture	29
5.1.2.2	Compensation of the floater with distilled water.....	29
5.1.3	Printing the value of ASV and of relative density	30
5.1.4	Zeroing of the balance	31
5.2	TOTAL DRY EXTRACT	31
5.2.1	Printing the value of TOTAL DRY EXTRACT	33
5.3	EVALUATION OF MUSTS	34
5.3.1	Printing the value of the degrees Baumé, Babo, Brix and Oechsle of a must.....	35
6.	SETTINGS OF THE SUPER ALCOMAT	36
6.1	Modification of the progressive label.....	36
6.2	Adjusting the date	36
6.3	Adjusting the time	36
6.4	Setting the contrast of the Super Alcomat display	36
6.5	Language.....	37
6.6	SERVICES (SUPPLEMENTARY CALIBRATION)	37
7.	SERIAL TRANSMISSION.....	38
7.1	CONNECTION OF THE SERIAL OUTPUT	38
7.2	RS232 COMMUNICATION PROTOCOL BETWEEN THE SUPER ALCOMAT AND THE PC	38
8.	USE AND MAINTENANCE SUGGESTIONS.....	39
9.	SIT CALIBRATION	39

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1. TECHNICAL CHARACTERISTICS SUPER ALCOMAT

GENERAL CHARACTERISTICS

Power tension	100/240 Vac -15/+10% (85÷264V) 50/60Hz
Power consumption	10 VA
Operating temperature	10 ÷ 30 °C (suggested 15 ÷ 25 °C)
Dimensions	210 x 370 x 380 mm (L x D x H)
Net weight	8,5 kg
Data output	RS232 I/O
Standard equipment	Power pack 2 hydrorepellent floaters, interchangeable in mass and volume Container for exhaust liquids - Certified hydroalcoholic solution

HYDROSTATIC BALANCE

Measurement range for the relative density	0,5 ÷ 2,25
Sensibility	0,00005
Reading repeatability	± 0,00005
Autocalibration with internal mass	YES
Response time	6 sec.

THERMOMETER

Type	Thermo probe in platinum Pt100 1/3 DIN
Maximum temperature displayed	50 °C
Sensibility	0,05 °C
Reading repeatability	± 0,05 °C (in the range 10 ÷ 30 °C)
Readability	0,1 °C

ALCOHOLIC STRENGTH BY VOLUME (ASV)

Reading range	0,13 ÷ 99,98 %ASV
sensibility	0,01 %ASV
Accuracy and reading repeatability	± 0,03 %ASV (between 15 and 25 °C)

TOTAL DRY EXTRACT

Range of values	0,0 ÷ 505,8 g/l of sucrose
Accuracy and reading repeatability	± 0,3 g/l

MUSTS EVALUATION

Relative density $d_{20/20}$ range	1,04 ÷ 1,16
------------------------------------	-------------

ACCESSORIES ON REQUEST:

- Floaters with certificate
- SIT certified solutions 5 ÷ 60 %ASV

FLOATERS INTERCHANGEABLE IN MASS AND VOLUME

A particular care is devoted to the construction and calibration of the floaters.

They are all interchangeable in mass and volume to simplify their use; for example, it is possible to use two different floaters to determine the relative density of a must or a wine; in this way possible sample pollution is avoided.

In order to achieve what mentioned before, the suspension wire¹ and the superficial treatment of the floater's glass were object of thorough studies to nullify the problems caused by the surface tension, by the wire's own volume, and by the electrostatic charges.

The floaters volume can, on request, be certified by competent organisations, such as the L.N.E. Paris, the I.M.G.C. Turin, the PTB Germany, etc.

¹ realised with a non-magnetic and stainless material, with a diameter less or equal to 0.20 mm.

2. INSTALLATION OF SUPER ALCOMAT

Super Alcomat was projected and built to resist and operate also in harsh conditions. However it is advisable to install the Super Alcomat:

- on a rigid base without vibrations
- in an environment with stable temperature and without excessive ventilation
- using a dedicated electrical line

2.1 INSTALLATION PROCEDURE

- Level the balance acting on the front supports so that the bubble is in the centre of the circle on the level;
- unlock the balance by rotating clockwise the knob on the right side of the instrument until it stops (Fig. 1 e Fig. 3 nr.15);


Fig. 1

- hang the floater on the dedicated prop (Figure 3 nr. 7);
- connect the circular plug of the transformer to the connector on the Super Alcomat indicated in Fig. 2 with nr.4. Insert then the plug in the network socket;

- 1 - Level bubble
- 2 - Printer connector
- 3 - RS232 connector
- 4 - Power socket
- 5 - Serial number
- 6 - Service connector


Fig. 2: Rear view

3. DESCRIPTION OF SUPER ALCOMAT

The instrument is constituted by a hydrostatic balance, with sensibility to the fifth decimal numeral, and an elaborating module integrated in the balance itself.



Fig. 3

- | | |
|-----------------------------------|---|
| 7- Floater prop | 13- Knob of the cylinder centering device |
| 8- Drainage pipe | 14- Cylinder centering device |
| 9- Cylinder and floater | 15- Knob for blocking balance |
| 10- Electrovalve | (to transport) |
| 11- Command panel/Display | 16- Thermo probe |
| 12- Supports for level regulation | 17- Tank for cylinder filling |
| | (minimum 65 ml) |

3.1 COMMAND PANEL

The keyboard of the Super Alcomat (Fig. 3 nr. 11) is equipped with a display and eight functional keys, to which a different meaning, depending on the displayed menu, is associated.

Every key must be pressed until an acoustic signal is uttered. This indicates that the instrument accepted the command.

NOTE

After some minutes when the Super Alcomat is not used, the display shuts off automatically. By pressing [ON], the display turns on again.

The keys allow the execution of the displayed action. The [MEASURE] and [DRAIN] keys activate always dedicated functions.

The [MEASURE] key allows:

- to automatically lower the floater for the determination of the specific gravity and of the ASV;
- to lift and lower the floater (press 1 sec.);
- to lift the floater the necessary time to allow an automatic zero setting of the system (press 3 sec).

The [DRAIN] key allows:

- to lift and hold in position the floater (press 1 sec.);
- to lift the floater and, at the same time, to open an electrovalve that allows the automatic emptying of the cylinder (press 3 sec). A flashing led and an acoustic signal indicate the opening time of the valve. After the emptying, the valve closes automatically.

4. IGNITION OF SUPER ALCOMAT

At ignition of the Super Alcomat, the software release number is displayed (Figure 3), anticipated by an acoustic signal and by the name of the company.


Fig. 3

After about two seconds, the MAIN MENU appears.


Fig. 4

By pressing the function keys [4] or [5] it is possible to access to "Settings" or "Services" (see paragraph 6); by pressing the function keys [1] or [2] or [3], the display shows the following blinking message:


Fig. 5

Any key pressed in this situation (Fig.5) brings back to the main menu, i.e. **the Super Alcomat during thermal stabilization cannot be used to perform measurements**, but the keyboard remains active.

The blinking message, if the keyboard is not touched, remains on the display until the Super Alcomat is in the thermal stabilization phase (about 5 minutes).

At the end of the thermal stabilization phase, the following blinking message appears on the Super Alcomat display:

**SELF CALIBRATION ACTIVATED
(PHASE 1) PLS. WAIT**

followed by

**SELF CALIBRATION ACTIVATED
(PHASE 2) PLS. WAIT**

The display of the Super Alcomat shows the MAIN MENU.

The system is now ready to work.

ATTENTION!

If at ignition the Super Alcomat displays the blinking message "OUT OF RANGE" press the key that indicates [-] on the Super Alcomat to continue. The blinking message "SCALE PERFORMING THERMAL STABILIZATION, PLS. WAIT" is displayed. At the end of the thermal stabilization, the Super Alcomat is operating. In case the thermal stabilization message persists for more than 15 minutes, press again [-] and have the Super Alcomat checked by a technician.

5. MEASUREMENTS WITH SUPER ALCOMAT

When the Super Alcomat is ready to work, the display shows the MAIN MENU (see paragraph 4, Figure 4). The MAIN MENU presents now the different functions that can be performed by the Super Alcomat.

By pressing one of the keys located on the inferior part or on the right side of the display, to which [1] o [2] o [3] o [4] o [5] correspond on the display, **the access to the desired corresponding function is granted**. The single functions are then guided by the indications on the display.

NOTE

For ease of writing, in the following text, the generic functional key correspondent to a certain text is indicated by the text displayed: for example, instead of saying "press the key correspondent to [EXIT]", it will be said "press [EXIT]".


Super Alcomat can operate at a temperature also **different than 20° (as long as in the range 15-25 °C)** and allows:

- through the function [1]: ALCOHOL. STRENGTH BY VOLUME (Alcoholic strength by volume or ASV)
 - to perform the relative density $d_{20/20}$ measurements of a distillate;
 - to visualise in 6 seconds and completely automatically the Alcoholic Strength by Volume of a distillate, from 0.13 to 99.98%ASV, in agreement with the official tables (EEC Regulations 2676/90 and 2870/00 which refer to the “International Alcoholometric Tables”, published in the International recommendation Nr. 22 of the O.I.M.L in 1972).
- through the function [2]: TOTAL DRY EXTRACT
 - to determine the total dry extract of a wine or of a must for specific gravity in the range 0.99 - 1.16;
- through the function [3]: MUSTS VALUATION
 - to evaluate the musts, by calculating the values of the Baumé, Brix and Oechsle degrees for $d_{20/20}$ in the range 1.04 - 1.16;

5.1 RELATIVE DENSITY AT 20°C OF A DISTILLATE AND ASV (ALCOHOLIC STRENGTH BY VOLUME)

It is done as follows:


- a- verify that the floater is in position on the dedicated prop;
- b- fill the cylinder with the liquid (minimum amount 65 ml), of which the specific gravity must be measured, by pouring in the apposite tank; the optimal replenishment level is granted by the automatic drainage systems which prevent the leak of the liquid itself; pay attention that no bubbles are present on the floater and in the liquid;
- c- press [1] to select “Alcoholic Strength by Volume”. The Super Alcomat display shows:

It indicates that the Super Alcomat is in the “zero condition”; it appears at the end of the calibration and/or in stand by condition. If anticipated by a “+” or a “-” it indicates that the hydrostatic balance is close to the zero.

Fig. 6: Zero condition of the hydrostatic balance

- d- press [MEASURE]: the floater will be automatically lowered; read the ASV value on the display
 - the ASV and the relative density at 20°C appear on the display after the measured relative density; the numbers change until the data does not reach stability;

Fig. 7: ASV and $d_{20/20}$ of a distillate

- a short acoustic signal and the ASV value shown with white fonts on black background highlight the achievement of this condition. The value remains on the display until the removal of the floater from the suspension;
 - the specific gravity shown is then related to 20°C in accordance with the tables required by the official method (EEC Regulations nr. 2676/90 and nr. 2870/00). The correction is reliable only for hydroalcoholic mixture and in the liquid temperature range 15 - 25°C. Outside this range, the results can undergo variations;
- e- if the analysis operations are completed, press [DRAIN] to empty the cylinder.

5.1.1 Relative density of a liquid

When the specific gravity of the liquid in analysis does not fall within the range of values that correspond to an Alcoholic Strength by Volume, the display shows only the relative density of the liquid and the temperature (Figure 8):


Fig. 8: Relative density of a liquid

5.1.2 Correction of the floater or compensation

Through the procedure of "Floater compensation" both the mass and the volume of the floater in use **are not modified**.

The so-called "compensation factor" of the floater is modified instead, so that all the floaters used will give the same response in the same standard conditions. These conditions are the determination of the relative density of distilled water and the determination of the ASV or alcoholic strength of a hydroalcoholic mixture, of known ASV and SIT certified, provided on request by our company. It is possible to select the two possibilities with the keyboard, by pressing [H₂O] or [ALCOHOL] respectively.

The display shows the current value of the compensation factor that can be increased or decreased by five units at a time by pressing [+] and [-]; the compensation factor remains stored even if the instrument is turned off.

It is advisable to check regularly this factor.

5.1.2.1 Compensation of the floater with the Gibertini SIT certified hydroalcoholic mixture

It is done as follows:

- verify that the floater is in position on the dedicated prop;
- fill the cylinder with a hydroalcoholic mixture of known strength (minimum amount 65 ml), at a temperature between 15 and 25°C; pay attention that no bubbles are present on the floater and in the liquid;
- select "Alcoholic Strength by Volume" from the main menu by pressing [1] (Figure 4), access to the option "Compensation" with the corresponding key (Figure 7);
- the display shows the current compensation factor, the relative density of the hydroalcoholic mixture, the specific gravity of the same mixture corrected at 20°C and the alcoholic strength (Figure 9);
- if the value of the alcoholic strength shown by the display is equal to that declared on the certificate of the mixture used as reference, then the instrument is in the ideal reading conditions; on the contrary, it is necessary to vary the compensation factor of the floater. Press [+] or [-], keeping in mind that [+] increases the alcoholic strength read, and [-] decreases it. It is necessary to point out that the variation of the compensation factor and the variation of the alcoholic strength don't correspond exactly. Moreover these variations are very slow: proceed therefore waiting some seconds before reading the value on the display, and stop on the value that is closest to the reference value.

The ASV measured with the Super Alcomat and that on the mixture certificate must agree


Decrease the ASV

Increase the ASV

ASV or alcoholic strength calculated starting from the relative specific gravity read

Relative density of the hydro-alcoholic mixture read with Super Alcomat

Previous screen

"Compensation of the floater in water"

It allows a manual zeroing, before the reading

Fig. 9: Compensation of the floater with the hydroalcoholic solution

5.1.2.2 Compensation of the floater with distilled water

It is done as follows:

- verify that the floater is in position on the dedicated prop;
- fill the cylinder with distilled water (minimum amount 65 ml) at a temperature between 15 and 25°C; pay attention that no bubbles are present on the floater and in the liquid;
- select "Alcoholic Strength by Volume" on the main menu by pressing [1] (Figure 4), then press [COMP.] to access the "Compensation" option (Figure 7) and press [H2O] (Figure 9);
- the display shows the current compensation factor, the relative density that water should have at the measured temperature, and the relative density of the water actually measured by the instrument (Figure 10);

- e- if the two densities agree, the instrument is in the ideal reading conditions; on the contrary, it is necessary to vary the compensation factor of the floater. Press [+] or [-], keeping in mind that [+] increases the alcoholic strength read, and [-] decreases it.


Fig. 10: Compensation of the floater with distilled water

5.1.3 Printing the value of ASV and of relative density

To have a printout of the value of the ASV or of the relative density of the liquid under consideration, just press [PRINT]; the option is available on the display at the moment of the determination. If the printer is not connected, the display shows the following blinking message:

**PRINTING LABEL PLEASES WAIT
PRINTER ERROR
PRESS ANY KEY TO CONTINUE**

If the printer is activated before a stable value of dry extract is reached, the display shows the blinking message:

UNSTABLE DATA: OPERAT. IMPOSSIBLE

While printing, the display of the Super Alcomat shows the blinking message:

**PRINTING LABEL PLEASES WAIT
PRINTING OK**

The printout is so composed, in the two possible cases:

Date: 06-03-2001 Hour : 16:07 Pr. 012
Alcoholic strength: [%ASV] : 13.24

Temp. of distillate : 23.7 °C
Spec. gravity at 20 °C : 0.98270
Spec. gravity of distillate : 0.98170

Date: 06-03-2001 Hour : 15:20 Pr. 001

Temp. of liquid : 24.5 °C
Spec. gravity of liquid : 1.10510

NOTE

On the label printed an identification number appears, indicated with "Pr", i.e. "progressive label". When the instrument is turned on, the numbering begins with 1. Every time a printout is made, the number is updated. It is possible to modify the numbering. For details, see paragraph 6.1 Modification of label number.

5.1.4 Zeroing of the balance

The instrument offers the possibility to perform an automatic zeroing during the measuring procedure. By keeping pressed [MEASURE], the floater is lifted for some seconds, allowing the balance to perform the zeroing. Super Alcomat offers also the possibility to perform a manual zeroing of the balance:

- a- if the cylinder is full, press [DRAIN], to drain the sample and lift the floater;
- b- press [ZERO]; the display shows “zero” followed by the initial screen;
- c- if the suspension is occupied, the display shows:

ZEROSET IMPOSSIBLE

5.2 TOTAL DRY EXTRACT

The total dry extract includes all matter that is non-volatile under specified physical conditions. The official European method is densimetric (EEC Regulation Nr. 2676/90 Attachment 4: Total Dry Extract). The total dry extract is calculated indirectly basing upon the value of the specific gravity of the alcohol-free must/wine, i.e. of the must/wine from which the alcohol was eliminated, and which was brought back to the initial volume by adding distilled water.

The total dry extract is expressed by the concentration in g/l of a sucrose solution with the same specific gravity of the alcohol-free must/wine.

It is done as follows:

- a- verify that the floater is in position on the dedicated prop;
- b- select TOTAL DRY EXTRACT from the page MAIN MENU, by pressing [2] (Figure 4). The instrument awaits for the measure of the relative density of the wine/must as it is (valid samples have specific gravity between 0.99 and 1.16) with the blinking message:


Fig. 11

- c- wait for the “zero” of the balance, or press [ZERO] for manual zeroing and press [MEASURE];
- d- the display shows temperature, specific gravity read, specific gravity correct at 20°C and density at 20°C of the wine/must sample under consideration (Figure 12);


Fig. 12

Super Alcomat relates automatically the specific gravity value read to 20°C.

Once a stable value is reached, the specific gravity value correct at 20°C (expressed with four decimal numerals) remains fixed on the display, with white fonts on black background, preceded by an acoustic signal.

The density is obtained by multiplying the specific gravity value correct at 20°C for 0,998203 g/ml (density of water at 20°C);

IT IS NOW POSSIBLE TO STORE UP TO 99 SAMPLES OF WINE/MUST

This function allows performing all the measurements on the samples of wine/must *in advance*, to recall them afterwards, at the moment of measuring the specific gravity of the correspondent distillates. To prevent sample pollution, use another floater (given that the compensation factor is correctly set, see paragraph 5.1.2);

To store a sample:

- press [SAVE] (Figure 12);
- through the keys [+], [-], [+10], [-10] assign to every sample an identification number;
- confirm the value with [SAVE] (operation allowed only if the value of relative density at 20°C is stable) and pass to the measurement of another sample of wine/must by emptying the cylinder.
- drain the previous liquid by pressing [DRAIN];
- wash abundantly with distilled water, in order to accurately clean the cylinder from any remainder. Rinsing with wine is suggested to prevent contamination, in case a sufficient amount of liquid is available; pour the alcoholic distillate of wine/must, wait for the "000" of the balance, press [MEASURE] to proceed to the next measurement.

To recall a sample:

- press [RECALL];
- through the keys [+], [-], [+10], [-10] select the identification number of the sample to recall;
- confirm the selection by pressing [RECALL] and pass to the determination of the dry extract by pressing [CALCULAT].

- e- to pass directly to the determination of the dry extract, press [CALCULAT]. The instrument awaits to measure the relative density of the distillate obtained by the wine/must under consideration (valid samples have a relative density $d_{20/20}$ between 0.79000 and 0.99995), with the blinking message (Fig. 13):


Fig. 13

- f- press [MEASURE]; the display shows temperature, specific gravity read and Alcoholic Strength by Volume of the distillate, together with the value of the total dry extract in g/l of sucrose (Figure 14);


Fig. 14

Super Alcomat provides automatically to determine the ASV of the distillate, using the correction algorithm valid for hydroalcoholic mixtures, and successively the relative density $d_{20/20}$ of the distillate itself.

Super Alcomat applies Tabarié's formula to determine the relative density $d_{20/20}$ of the alcohol free wine/must, then refers the correspondent tabulated value⁵ of total dry extract.

Once a stable value is reached, the total dry extract value is fixed on the display, with white fonts on black background, preceded by an acoustic signal;

- g- once the total dry extract value is calculated, it is possible to print the relative label (see next paragraph) and/or delete from the memory the sample of wine/must just analysed by keeping pressed [DELETE];
- h- to exit the total dry extract function press [EXIT] until the display shows the MAIN MENU;
- i- at the end of the measurement, drain the liquid analysed by pressing [DRAIN];
- j- at the end of every measurement, wash abundantly with distilled water, in order to accurately clean the cylinder from any remainder.

5.2.1 Printing the value of TOTAL DRY EXTRACT

To print the value of TOTAL DRY EXTRACT of the wine/must under consideration, just press [PRINT]; the option is available on the display at the moment of the determination.

While printing, the display of the Super Alcomat shows the blinking message:

PRINTING LABEL PLEASES WAIT
PRINTING OK

⁵ For further details, see R. Becchetti, "Analysis methods of wines and alcoholic mixture", 1999 by Gibertini elettronica

The label is so composed:

Date: 09-03-2001 Hour : 10:24	
Sample N. : 29●	

Dry extract [g/l]	: 129.2
Alcoholic str. [%ASV]	: 3.95
Spec. Gravity wine/must at 20°C:	1.044

The identification number of the sample is chosen during the storing phase (see previous paragraph)
 "Sample Nr." IS NOT updated of 1 unit after every printout and it is the last value selected in [SAVE].

If the printer is not connected, the display shows the following blinking message:

**PRINTING LABEL PLEASES WAIT
 PRINTER ERROR
 PRESS ANY KEY TO CONTINUE**

If the printer is activated before a stable value of dry extract is reached, the display shows the blinking message:

UNSTABLE DATA: OPERAT. IMPOSSIBLE

5.3 EVALUATION OF MUSTS

The degrees Baumé, Babo, Brix and Oechsle are in relation with the relative density $d_{20/20}$ of a must. The use of the Super Alcomat allows to determine automatically the density and the value of the degrees Baumé, Babo, Brix and Oechsle of a must, for density $d_{20/20}$ between 1.04 and 1.16, in the range of temperature 15°C - 25 °C. Below 15°C and above 25 °C the results are not completely reliable.

It is done as follows:

- in case the must is not clear, perform a filtration;
- verify that the floater is in position on the dedicated prop;
- select MUSTS EVALUATION in the MAIN MENU screen, by pressing [3] (Figure 4). Wait for the "zero" of the balance and press [MEASURE]. The floater is automatically submerged in the liquid and the measurement begins. The display shows temperature, specific gravity read, specific gravity correct at 20°C (with four decimal numeral) and density at 20°C of the must sample under consideration (the sample is valid for $d_{20/20}$ between 1.04 and 1.16, Figure 15);
- after some seconds, after a stable relative density value is reached, the numerical values correspondent to the degrees Baumé, Babo, Brix and Oechsle of the must are shown on the display, preceded by a short acoustic signal.


Fig. 15

Super Alcomat determines the relative density $d_{20/20}$ of the must under consideration. The density of the must is obtained by multiplying the value of the specific gravity correct at 20°C for 0,998203 g/ml (density of water at 20°C).

In correspondence of the density correct at 20 °C, the Super Alcomat⁶

- reports the degree Babo, Brix, etc. correspondent if the values of density are present in the table;
- performs an interpolation for values of density not present in the table;

IT IS POSSIBLE TO ENUMERATE UP TO 99 MUST SAMPLES

- press [NUMBER];
- by pressing [+], [-], [+10], [-10] assign to every sample an identification number and press [EXIT];
- the chosen numerical value is shown on the display and on the possible printout. It is automatically updated every time a printout is made. It is possible therefore to identify the sample under consideration through the function [NUMBER] or through the automatic updating of 1 unit of the sample enumeration.

e- at the end of the measurements, drain the liquid by pressing [DRAIN];

f- wash abundantly with distilled water, in order to accurately clean the cylinder from any remainder;

g- press [EXIT] to go back to the MAIN MENU.

5.3.1 Printing the value of the degrees Baumé, Babo, Brix and Oechsle of a must

To have a printout relative to the evaluation of a must, just press [PRINT]; the option is available on the display at the moment of the determination.

While printing, the display of the Super Alcomat shows the blinking message:

**PRINTING LABEL PLEASES WAIT
PRINTING OK**

The label is so composed:

Date: 10-03-2001	Hour	:10:38
Musts evaluation		
Must sample N.	:	16●
Temperature of must	:	24.6 °C
Spec. gravity of must at 20 °C	:	1.1060
Degr. Baumé	:	13.93
Degr. Babo	:	21.35
Degr. Brix	:	24.9
Degr. Oechsle	:	106.0

Number automatically increased of 1 unit every time a printout is performed or value chosen with [NUMBER].

If the printer is not connected, the display shows the following blinking message:

**PRINTING LABEL PLEASES WAIT
PRINTER ERROR
PRESS ANY KEY TO CONTINUE**

If the printer is activated before a stable value of density is reached, the display shows the blinking message:

UNSTABLE DATA: OPERAT. IMPOSSIBLE

⁶ For further details, see R. Becchetti, "Analysis methods of wines and alcoholic mixture", 1999 by Gibertini elettronica

6. SETTINGS OF THE SUPER ALCOMAT

To modify some settings, the options menu presents different possibilities that are described in the following.

6.1 Modification of the progressive label

When determining the ASV of a distillate, or the relative density of a liquid, it is possible to identify the sample under consideration through a number to report while printing. This number is represented by the progressive label, which, as already said, increases of 1 unit after every printout, and the ASV or density value is different from the previous one.

To modify the "Progressive label" number, select from the MAIN MENU "Settings", by pressing [4] (Figure 4), then select "Progressive label" by pressing [2]. By pressing [+] and [-] that progressive number is varied of 1 unit, while with [+10], [+100] tens and hundreds are selected. The maximum possible value for the labels is 999.

The labels that include only the value of the specific gravity and the temperature have a numbering that cannot be varied by the user.

6.2 Adjusting the date

To adjust the date, select "Settings" from the MAIN MENU by pressing [4] (Figure 4), then select "Set the date" by pressing [3].

Highlight the day or the month or the year to modify by pressing respectively [DAY] or [MONTH] or [YEAR]. At this point modify the value by pressing [+] or [-] keeping in mind that the two keys increase and decrease respectively of 1 unit the value under adjustment.

6.3 Adjusting the time

To adjust the time, select "Settings" from the MAIN MENU by pressing [4] (Figure 4), then select "Set the time" by pressing [4].

Highlight the hour or the minutes or the seconds to modify by pressing respectively [HOUR] or [MIN] or [SEC]. At this point modify the value by pressing [+] or [-] keeping in mind that the two keys increase and decrease respectively of 1 unit the value under adjustment.

6.4 Setting the contrast of the Super Alcomat display

To adjust the contrast of the Super Alcomat display, select "Settings" from the MAIN MENU by pressing [4] (Figure 4), then select "Display Contrast" by pressing [5].

The [REVERSE] key (Figure 16) allows inverting the contrast: writings become white on black background, while [NORMAL] brings back to the initial contrast: black writings on white background. By pressing [+] the display background becomes darker, while [-] makes it lighter.


Fig. 16

6.5 Language

To change the language used to express all the messages that appear on the display, select “Settings” from the MAIN MENU by pressing [4] (Figure 4), then select “Language” by pressing [6]. Identify the desired language (Figure 17) among Italian, English, French, Spanish and German by pressing the correspondent key. Press [1] to exit.


Fig. 17

All the messages that appear on the display will be shown in the new language chosen. The selection will be maintained also if the instrument is turned off.

6.6 SERVICES (SUPPLEMENTARY CALIBRATION)

It must be remembered that the instrument provides automatically, after about 10 minutes at rest, to start the calibration procedure.

If the operator senses it necessary, he can activate a supplementary calibration.

Select “Services” from the MAIN MENU by pressing [5] (Figure 4), then select “Calibration” by pressing [2].

The display of the Super Alcomat shows:

**CALIBRATION PROCEDURE ACTIVATED
PLS. WAIT**

then the blinking message:

**SELF CALIBRATION ACTIVATED
(PHASE 1) PLS. WAIT**

followed by

**SELF CALIBRATION ACTIVATED
(PHASE 2) PLS. WAIT**

During the first and the second calibration phases it is not possible to use the instrument.

7. SERIAL TRANSMISSION

The serial transmission of the SUPER ALCOMAT uses the standard EIA RS232.

7.1 CONNECTION OF THE SERIAL OUTPUT

The connection is the following:

Signals	Connector DB9F SUPER ALCOMAT	Connector DB9F HOST
RS 232 Out	Pin 2	Pin 3
RS 232 In	Pin 3	Pin 2
Ground	Pin 5	Pin 5

7.2 RS232 COMMUNICATION PROTOCOL BETWEEN THE SUPER ALCOMAT AND THE PC

By means of the RS232 communication protocol, it is possible to receive with a computer the data visualised on the Super Alcomat according to the following protocol:

- Configuration of the serial port COM1 : 9600, 8, N, 1.

COMMANDS TABLE

Valid Commands	Answer to "legal" commands	Answer to "illegal"* commands
T<CR>	<Temperature><CR>	<----><CR>
S<CR>	<Density read><CR>	<----><CR>
D<CR>	<Density corrected at 20 °C><CR>	<----><CR>
A<CR>	<Alcoholic strength><CR>	<----><CR>
E<CR>	Phase 1: <Dry extract>, <----><CR> Phase 2: <Dry Extract>, <Alcoholic Strength><CR>	<----><CR>
M<CR>	<Mass>, <Baume>, <Babo>, <Brix>, <Oe><CR>	<----><CR>

*- The command is defined as illegal when the operating condition does not allow its visualisation.

N.B.- <CR> means the ASCII character "Carriage Return"

SUPER ALCOMAT provides always a correct answer to all the commands received

- if the command is recognised, and the context is valid, the answer will be the value requested, followed by the end-line character **<CR>**;
- if the command is correctly recognised, but it is sent out of a valid context (for example: request of the total dry extract while measuring alcohol), the answer will be a series of 4 dashes **<---->** followed by the end-line character **<CR>**;
- if the command is not recognised (different command than those indicated in the table) the answer will be a **< ? >** followed by the end-line character **<CR>**.

8. USE AND MAINTENANCE SUGGESTIONS

It is advisable to group the hydroalcoholic solutions to examine in homogeneous groups, for example: from 9% to 12% ASV - from 15% to 20% ASV - from 40% to 45% ASV.

This is made to prevent pollution between a determination and the following.

In other words, it is **NOT** advisable to pass from 10% to 40% and vice versa, unless after doing an abundant rinsing with wine, if enough liquid is available. When passing from 10% to 40%, wash profusely with H₂O to prevent "contamination".

Wash with H₂O or alcohol at the end of the work.

A dedicated microprocessor provides, besides controlling the most important functions of the Super Alcomat, to diagnose internal problems of the balance, which may have effects on the precision, suggesting a programmed service after 10000 hours of functioning with qualified personnel. The displayed message is "**The instrument requires maintenance press a key to continue**" which, however **does not** compromise the use of the balance.

9. SIT CALIBRATION

Gibertini Elettronica is SIT CENTER n. 94 (law 273/91) accredited for the characterization of electronic balances with capacity till 30 kg.

The certifications issued have an official value: they are documents that, according to the law, can be used for all the procedures of certification, homologations and accreditation of products, services and systems of business quality.

Gibertini, therefore, can emit the calibration certificate that guarantees:

- the maintenance of the referral of apparatus used by the Centre with national samples of unity of the International System of Units (S.I.);
- the metrological correctness of the measuring procedure adopted by the Centre.

UPON REQUEST: SIT Certified hydroalcoholic solutions from 5%ASV to 60%ASV.

MAPS OF WORLDWIDE RECOGNIZED INTERNATIONAL CENTRES EQUIVALENT TO SIT

NOTE

info@gibertini.com
Tel. (++39) 023541434

ISO 9001:2008

LAT N° 094

**GIBERTINI
ELETTRONICA s.r.l.**

Via Bellini 37
20026 Novate Mil.se
Milano (Italy)

Tel. (++39) 02 3541434
Fax (++39) 02 3541438
www.gibertini.com
sales@gibertini.com